

Instruction



Bend this lip
over (min 1 mm)





Reason: The DFI panel is bent under pressure of this lip.

Tools: Combination Pliers